

Approved by Director: Dr. Guy Vallaro

DAY 1

Wear a clean disposable lab coat

Clean bone pulverizing equipment:

Locate cylinder, end pieces, and impactor

Scrub equipment using 20% (in-house) bleach or 10% stabilized bleach, warm Tergazyme, and/or DNA-Off

Place equipment into 50mL Falcon tubes containing DNA-Off or bleach for 20 minutes on the rotator

Rinse parts thoroughly with dH₂O

After rinsing, fill Falcon tubes with dH₂O and place back on rotator for 20 minutes

Dry down equipment using isopropyl wipes

Stratalink equipment for at least 15 minutes (45 minutes is optimal)

Clean dead space hood:

Fresh 20% (in-house) bleach or 10% stabilized bleach

Isopropyl wipe 1 pipettor (p1000)

Parafilm strips

UV for at least 15 minutes

Clean Air Clean Hood (bone/tooth cutting):

Fresh 20% (in-house) bleach or 10% stabilized bleach

Bleach/isopropyl wipe Dremel

UV for at least 15 minutes

Pipette:

1.7mL EDTA (850 μ L x2, kept over sink) into labeled extraction tube(s)

Clean:

Dremel bit, forceps, and scissors (if used) using bleach and isopropyl wipes

Stratalink for at least 15 minutes:

6 extraction tubes (2.0mL) labeled RB, Q, and 4 powder tubes

*adjust the number of labeled tubes if more than one RB tube will be extracted

Extraction tube(s) (2.0mL) with 1.7mL EDTA

15mL Falcon tube for RB's

2 cotton swabs

Approved by Director: Dr. Guy Vallaro

Cutting disc
2 sanding discs
Dremel bit
Screwdriver
Scissors (if used)
Forceps
2 sheets of weigh paper folded diagonally

In dead space hood:

Move everything from the Stratalinker to the dead space hood
Swab the inside of the cylinder, end pieces, and impactor
Cut ~1/2-1/3 of swab into RB tube, place remainder in 15mL Falcon tube (snap off shaft)
Add 1.6mL EDTA to RB swab (800µL x2)
Add 1.6mL EDTA to additional RB tubes as necessary (800µL x2)
Parafilm tube(s), place on room temperature nutator overnight
Repeat swabbing of parts with 2nd swab, place entire swab into the 15mL Falcon tube (snap off shaft)
Store tube of swabs at 4°C (short term) or -20°C (long term)
Assemble the cylinder, end pieces, and impactor

Transfer:

Dremel bit, screwdriver, and sanding/cutting discs to Air Clean Hood

Take out evidence:

Photograph evidence
Document a description on Worksheets

In Air Clean Hood:

Wear goggles and change gloves immediately after cutting

Place bone/tooth in hood
Screw sanding disc onto Dremel bit
Sand the exterior surface of the bone/tooth (use the second disc if necessary)
Screw the cutting disc on

For a tooth: Remove the root(s) (cut at the dentin)

For a bone: Remove ~1-2cm X 1-2cm section (make a wedge-shaped cutting)

*Remember to take a photograph of the evidence after you take a cutting- this includes any soft tissue and/or bone marrow that is removed from the evidence.

Approved by Director: Dr. Guy Vallaro

Pulverize sample:

Transfer cylinder, end pieces, and impactor (already assembled) to the counter
Place excised portion of sample into cylinder under the impactor
Place the cylinder into the freezer/mill

Grab the Cryogenic gloves, and keep on the goggles

(Liquid Nitrogen located downstairs in the Chemical room)

Fill the freezer/mill with liquid Nitrogen to the “fill line”
Close the freezer/mill and allow the liquid Nitrogen to freeze the sample for 10 minutes
After 10 minutes, open the freezer/mill and re-fill with liquid Nitrogen to the “fill line”

Using the touch LED screen, set the program as follows: Pre cooling period- 1 minute, Running period- 5 minutes, Cooling period- 2 minutes, Rate 15 cps, Cycle- 1

Press anywhere on the main screen
Press on the “Settings” tab
Move the bar along the time line to the desired settings
Once done, press “Return” tab
Press “Start” tab

When the freezer/mill is done with the cycle, remove the cylinder (use the cryogenic gloves)
Open the cylinder with the extractor and transfer the powder onto the weigh paper
Place enough powder for a single extraction into each tube (~1/8- 1/4 of a 2.0mL tube)
Store powder tubes at 4°C (short term) or -20°C (long term)

Cleanup:

Bleach/isopropyl wipe Dremel
Place discs, plastic cylinder, and weigh papers into biohazard waste
Scrub Dremel bit, impactor, and end pieces with warm Tergazyme and/or bleach
Repeat initial DNA-Off or bleach / dH₂O cleaning of end pieces and impactor
Change pre-filter in Air Clean Hood
Bleach down Air Clean Hood and counters
Place disposable lab coat, pre-filter, and gloves in biohazard box
Bleach and clean out dead space hood

Perform at **END OF DAY** to avoid over-decalcification:

Clean dead space hood using 20% (in-house) bleach or 10% stabilized bleach
Isopropyl wipe pipettor (p1000)
Pipette 1.7mL EDTA into labeled extraction tube(s) (2.0mL)
Stratalink tube(s) of EDTA and piece(s) of parafilm
UV hood for at least 15 minutes

Approved by Director: Dr. Guy Vallaro

Add 1.6mL EDTA to powder tube(s) (800 μ L x2)
Parafilm tube(s) and invert to thoroughly mix the liquid and powder
Place on room temperature nutator for 8-16hrs (up to 24hrs acceptable)

DAY 2- Morning

Clean Laminar Flow Hood: Fresh 20%(in-house) bleach or 10% stabilized bleach
Isopropyl wipe 3 pipettors (p2, p200, p1000)
UV at least 15 minutes

Make SEB/DTT: 5mL SEB (pre-aliquoted) + 0.03g DTT (refrigerator) in 15mL Falcon tube

Stratalink: SEB/DTT, dH₂O in 15mL Flacon tube, 1.5mL rack for at least 15 minutes

Extraction: Remove parafilm from RB(s), place in centrifuge
Repeat for Q(s) Spin for 1 minute at 8000 rcf

For the remainder of the procedure:

Manipulate RB tube(s) and place into centrifuge before touching Q tube(s)
Change gloves each time after handling Q tube(s)
Apply UV to Laminar Hood during spins

Pipette off EDTA from RB(s) using p1000/p200
Add 1mL dH₂O, place in centrifuge
Repeat for Q(s), invert/flick to re-suspend (do not vortex)
Repeat spin

Pipette off dH₂O from RB(s) using p1000/p200
Add 1mL dH₂O, place in centrifuge
Repeat for Q(s), invert/flick to re-suspend (do not vortex)
Repeat spin

Pipette off dH₂O from RB(s) using p1000/p200
Add 1mL dH₂O, place in centrifuge
Repeat for Q(s), invert/flick to re-suspend (do not vortex)
Repeat spin

Pipette off dH₂O from RB(s) using p1000, p200
Add 300 μ L SEB/DTT, 2 μ L proK to RB(s)
Repeat for Q(s), invert/flick to re-suspend
Place on nutator in 56°C incubator w/ FAN ON 6-8 hrs (up to overnight)

Approved by Director: Dr. Guy Vallaro

DAY 2- Afternoon (or DAY 3- Morning)

Stratalink: Rack with dH₂O

Microcon sets (filter and tube)

Additional Microcon tubes (according to how many RB and sample tubes you have).

Clean Laminar Flow Hood: Fresh 20%(in-house) bleach or 10% stabilized bleach

Falcon tube rack

Isopropyl wipe 2 pipettors (p200, p1000)

UV at least 15 minutes

For the remainder of the procedure manipulate the RB tube and place into centrifuge before touching Q tube. Change gloves each time after handling Q tube. Apply UV to hood during spins.

Pulse spin tubes that have been incubating to collect condensate

Add 300 µL PCIA to each tube, vortex, and spin for 3 minutes at 10,000 g (rcf)

Add 100 µL of dH₂O to Microcon set while waiting for tubes to spin.

Pipette off supernatant and add to Microcon set, spin for 5 minutes at 500-3,000 g (rcf)

Transfer filter to new Microcon tube, add 400 µL of dH₂O, and spin for 5 minutes at 500-3,000 g (rcf)

Add 60 µL of dH₂O to the filter, invert into a new Microcon tube, vortex, and spin for 3 minutes at 10,000 g (rcf).

Using a sterile pipette tip, determine the volumes of the RB and the sample extracts. The elution volumes shall be documented manually on QRM-5. The volume of the RB must not exceed the volume of the sample. If necessary, add dH₂O to bring the sample up to the volume of the RB.

Transfer RB extracts into one new Microcon set. Transfer the sample extracts into another Microcon set. Spin for 5 minutes at 500-3,000 g (rcf)

Add 60 µL of dH₂O to the filter, invert into a new Microcon tube, vortex, and spin for 3 minutes at 10,000 g (rcf).

Using a sterile pipette tip, determine the volumes of the RB and the sample extracts. The elution volumes shall be documented manually on QRM-5. The volume of the RB must not exceed the

Approved by Director: Dr. Guy Vallaro

volume of the sample. If necessary, add dH₂O to bring the sample up to the volume of the RB.

ARCHIVED